

U2 XML Example

Below illustrates some U2 XDOM functions for working with an XML document. The XML document should be placed in file &XML& with the ID of "SAMPLE.XML" to work with this code.

This example was created in response to the follow message posted on the u2-users group on 8th October 2008 which can be found at <http://listserver.u2ug.org/>.

```
Hello all,
```

I have to be able to parse out an XML file being sent by an laboratory instrument.

Below is a sample file... I need to be able to get to the ID attribute in SA, as well as the Key in AR, and AR's data.

I need to do this in BASIC, and pass the parsed data to another routine.

I've tried with OpenXMLData, ReadXMLData with an extraction file that I'm not positive how to build (samples I found don't deal with attributes at all).

I also tried with the XDOMOpen/Locate/etc. with no luck either. I'm about to go off and parse this thing myself but thought I'd try the list before I do.

Any help is appreciated!

TIA,

Robert

XML from the message above cleaned up

Save this as "SAMPLE.XML" in file &XML&

```
<?xml version="1.0" encoding="utf-8"?>
<?xml-stYLESHEET type="text/xsl" <?xml version="1.0" encoding="UTF-16"?>
<?xml-stYLESHEET type="text/xsl" href="C:\IRIS2K1\Templates\SpecimenAnalysis-style.xml"?>
<SA BF="URN" ID="021305941" SID="AP" OP="gloa" ADT="2005-06-08 16:03:32-08:00" ADTS="2005-06-08 16:03:32" RDT="2005-06-10 17:15:44-08:00" RDTs="2005-06-10 17:15:44" RP="2" SQN="0" RSQN="28" DILN="1" DILD="1" IMP="0" CDT="2005-06-08 16:03:32-08:00" CDTs="2005-06-08 16:03:32" REDT="2005-06-08 16:13:31-08:00" REDTS="2005-06-08 16:13:31">
  <PF></PF>
  <PF></PF>
  <PF></PF>
  <PF></PF>
  <PF></PF>
  <PF></PF>
  <PF></PF>
  <AC AT="Chemistry" AS="Done" SO="External">
    <AR Key="GLU" SN="GLU" LN="Glucose" AF="0" NR="30">Neg</AR>
    <AR Key="PRO" SN="PRO" LN="Protein" AF="0" NR="30">Neg</AR>
    <AR Key="BIL" SN="BIL" LN="Bilirubin" AF="0" NR="1+">Neg</AR>
    <AR Key="URO" SN="URO" LN="Urobilinogen" AF="1" NR="1+">3+</AR>
    <AR Key="PH" SN="PH" LN="pH" AF="1" NR="5.0">5.0</AR>
    <AR Key="BLD" SN="BLD" LN="Blood" AF="0" NR="Trace">Neg</AR>
    <AR Key="KET" SN="KET" LN="Ketone" AF="1" NR="1+">2+</AR>
    <AR Key="NIT" SN="NIT" LN="Nitrite" AF="1" NR="Pos">Pos</AR>
    <AR Key="LEU" SN="LEU" LN="Leukocytes" AF="1" NR="1+">4+</AR>
    <AR Key="CLA" SN="CLA" LN="Clarity" AF="0" NR="Hazy">Clear</AR>
    <AR Key="SG" SN="SPGR" LN="Specific Gravity" AF="0" NR="1.040">1.015</AR>
    <AR Key="COL" SN="COL" LN="Color" AF="0" NR="Amber">Colorless</AR>
  </AC>
  <AC AT="Sediment" AS="Done" SO="Internal">
    <AR Key="ART" SN="ART" LN="Artifact" AF="0" NR="99999999 /LPF">[none]</AR>
    <AR Key="RBC" SN="RBC" LN="Red Blood Cell" AF="1" NR="4 /uL">33 /uL</AR>
    <AR Key="WBC" SN="WBC" LN="White Blood Cell" AF="0" NR="6 /HPF">[none]</AR>
    <AR Key="WBCC" SN="WBCC" LN="White Blood Cell Clump" AF="1" NR="Occ">Many</AR>
    <AR Key="BACT" SN="BACT" LN="Bacteria" AF="0" NR="Few">Rare</AR>
    <AR Key="BYST" SN="BYST" LN="Budding Yeast" AF="0" NR="Few">[none]</AR>
    <AR Key="HYST" SN="HYST" LN="Hyphae Yeast" AF="0" NR="Few">[none]</AR>
    <AR Key="SQEP" SN="SQEP" LN="Squamous Epithelial" AF="0" NR="16 /HPF">[none]</AR>
    <AR Key="TREP" SN="TREP" LN="Transitional Epithelial" AF="0" NR="1 /HPF">[none]</AR>
    <AR Key="REEP" SN="REEP" LN="Renal Epithelial" AF="0" NR="1 /HPF">[none]</AR>
    <AR Key="OVFB" SN="OVFB" LN="Oval Fat Body" AF="0" NR="1 /LPF">1 /LPF</AR>
    <AR Key="FAT" SN="FAT" LN="Fat" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="MUCS" SN="MUCS" LN="Mucous" AF="0" NR="999999999 /LPF">[none]</AR>
    <AR Key="RBCC" SN="RBCC" LN="Red Blood Cell Clump" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="SPRM" SN="SPRM" LN="Sperm" AF="0" NR="[none]">Rare</AR>
    <AR Key="TRCH" SN="TRCH" LN="Trichomonas" AF="0" NR="Present">[none]</AR>
    <AR Key="NSE" SN="NSE" LN="Non-Squamous Epithelial" AF="0" NR="1 /HPF">&lt; 1 /HPF</AR>
    <AR Key="UNCC" SN="UNCC" LN="Unclassified Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="HYAL" SN="HYAL" LN="Hyaline Cast" AF="0" NR="3-5">[none]</AR>
    <AR Key="EPIC" SN="EPIC" LN="Epithelial Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="WBCT" SN="WBCT" LN="White Blood Cell Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="RBCT" SN="RBCT" LN="Red Blood Cell Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="GRAN" SN="GRAN" LN="Granular Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="CELL" SN="CELL" LN="Cellular Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="BROAD" SN="BROAD" LN="Broad Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="FATC" SN="FATC" LN="Fatty Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="WAXY" SN="WAXY" LN="Waxy Cast" AF="0" NR="1 /LPF">[none]</AR>
    <AR Key="UNCX" SN="UNCX" LN="Unclassified Crystal" AF="0" NR="1 /HPF">1 /HPF</AR>
    <AR Key="TPO4" SN="TPO4" LN="Triphosphate Crystal" AF="0" NR="FEW">[none]</AR>
    <AR Key="CAOX" SN="CAOX" LN="Calcium Oxalate Crystal" AF="0" NR="FEW">[none]</AR>
    <AR Key="CAPH" SN="CAPH" LN="Calcium Phosphate Crystal" AF="0" NR="FEW">[none]</AR>
    <AR Key="CACB" SN="CACB" LN="Calcium Carbonate Crystal" AF="0" NR="FEW">[none]</AR>
    <AR Key="URIC" SN="URIC" LN="Uric Acid Crystal" AF="0" NR="FEW">[none]</AR>
    <AR Key="LEUC" SN="LEUC" LN="Leucine Crystal" AF="0" NR="POS">[none]</AR>
    <AR Key="CYST" SN="CYST" LN="Cystine Crystal" AF="0" NR="POS">[none]</AR>
    <AR Key="TYRO" SN="TYRO" LN="Tyrosine Crystal" AF="0" NR="POS">[none]</AR>
    <AR Key="AMOR" SN="AMOR" LN="Amorphous Crystal" AF="0" NR="FEW">[none]</AR>
    <AR Key="UNCL" SN="UNCL" LN="Unclassified" AF="0" NR="99999999 /LPF">[none]</AR>
    <AR Key="PC" SN="PC" LN="PC" AF="0" NR="[none]">1314 /uL</AR>
  </AC>
  <FL>CHEMCONFIRM</FL>
  <CM>This sample is contaminated!! Comment appears here!!!</CM>
  <ARV>42</ARV>
</SA>
```

Example Code

```
$INCLUDE UNIVERSE.INCLUDE XML.H
*
    EQU TRUE TO 1
    EQU FALSE TO 0
*
    DIM    SA.AC.NODES(10)
    SA.NODE.CNT = 1
*
    RTN.CODE = XDOMOpen("SAMPLE.XML", XML.FROM.FILE, XDOM)
    IF RTN.CODE = XML.SUCCESS THEN
        CRT "XML Document opened"

        ;* read attribute "ID" from node "SA"
        XPATH = "/SA"
        ATT.ID = "ID"
        GOSUB 100
        IF NOT(ERR) THEN
            CRT "SA, Attribute (ID) = <":VALUE:>."
*
            ;*****
            ;* find all our nodes with a path of /SA/AC
            ;* and store them in our dimensioned array SA.AC.NODES
            ;*****
            XPATH = "/SA/AC"
            GOSUB 200
*
            IF NOT(ERR) THEN
                SA.AC.NODES(SA.NODE.CNT) = FND.NODE
                NEW.NODE = FND.NODE
                LOOP
                    RTN.CODE = XDOMLocateNode(NEW.NODE, XDOM.NEXT.SIBLING.WITH.SAME.NAME, 1,
XDOM.ELEMENT.NODE, NEW.NODE)
                    WHILE RTN.CODE = XML.SUCCESS
                        SA.NODE.CNT += 1
                        SA.AC.NODES(SA.NODE.CNT) = NEW.NODE
                        REPEAT
                    END
*
                    ;*****
                    ;* loop through our found nodes and find our "AR" children
                    ;*****
                    FOR X = 1 TO SA.NODE.CNT

                        ;* print out our AC, AT attribute value
                        FND.NODE = SA.AC.NODES(X)
                        ATT.ID = "AT"
                        GOSUB 300
                        CRT " " ":VALUE
*
                        ;* find our "AR" node
                        RTN.CODE = XDOMLocate(SA.AC.NODES(X), "AR", "", FND.NODE)
                        GOSUB 400
*
                        LOOP
                            RTN.CODE = XDOMLocateNode(FND.NODE, XDOM.NEXT.SIBLING.WITH.SAME.NAME, 1,
XDOM.ELEMENT.NODE, FND.NODE)
                            WHILE RTN.CODE = XML.SUCCESS
                                ;* locate our text node in "AR" to get results
                                GOSUB 400
                                REPEAT
                            NEXT X
                        END
                    END ELSE
                        CRT "ERROR: unable to open xml document!"
                    END
                RETURN
```

Example Code (cont)

```
*
100: *** read attribute from node ***
*
*   XPATH (IN):  path of the node to read
*   ATT.ID (IN): id of the attribute to get the value of
*
*   VALUE (OUT): value of the attribute, otherwise ""
*   ERR   (OUT): true if an error occurred
*
*   ERR = FALSE
*
*   VALUE = ""
*   GOSUB 200
*   IF NOT(ERR) THEN
*       GOSUB 300
*   END
*   RETURN

*
200: *** locate node ***
*
*   XPATH (IN):  path of the node to read
*
*   FND.NODE (OUT): node that was found
*   ERR      (OUT): true if an error occurred
*
*   ERR = FALSE
*   RTN.CODE = XDOMLocate(XDOM, XPATH, "", FND.NODE)
*   IF RTN.CODE # XML.SUCCESS THEN
*       CRT "ERROR: unable to locate XPATH <":XPATH:">"
*       ERR = TRUE
*   END
*   RETURN

*
300: *** read attribute ***
*
*   ATT.ID   (IN) : attribute to read
*   FND.NODE (IN) : node to read attribute from
*
*   VALUE    (OUT): value of the attribute, otherwise ""
*   ERR      (OUT): true if an error occurred
*
*   ERR = FALSE
*   RTN.CODE = XDOMGetAttribute(FND.NODE, ATT.ID, ATT.NODE)
*   IF RTN.CODE = XML.SUCCESS THEN
*       RTN.CODE = XDOMGetNodeValue(ATT.NODE, VALUE)
*   END ELSE
*       CRT "ERROR: unable to read attribute <":ATT.ID:">"
*       ERR = TRUE
*   END
*   RETURN

*
400: *** read results and display ***
*
*   ;* locate our text node in "AR" to get results
*   RTN.CODE = XDOMLocateNode(FND.NODE, XDOM.CHILD, XDOM.FIRST.CHILD, XDOM.TEXT.NODE, TEXT.NODE)
*   RTN.CODE = XDOMGetNodeValue(TEXT.NODE, TEXT.VALUE)
*
*   ATT.ID = "Key"
*   GOSUB 300
*   CRT "      ":VALUE:" = ":TEXT.VALUE
*   RETURN
*
END
```

Output

XML Document opened
SA, Attribute (ID) = <021305941>.

Chemistry

GLU = Neg
PRO = Neg
BIL = Neg
URO = 3+
PH = 5.0
BLD = Neg
KET = 2+
NIT = Pos
LEU = 4+
CLA = Clear
SG = 1.015
COL = Colorless

Sediment

ART = [none]
RBC = 33 /uL
WBC = [none]
WBCC = Many
BACT = Rare
BYST = [none]
HYST = [none]
SQEP = [none]
TREP = [none]
REEP = [none]
OVFB = 1 /LPF
FAT = [none]
MUCS = [none]
RBCC = [none]
SPRM = Rare
TRCH = [none]
NSE = < 1 /HPF
UNCC = [none]
HYAL = [none]
EPIC = [none]
WBCT = [none]
RBCT = [none]
GRAN = [none]
CELL = [none]
BROAD = [none]
FATC = [none]
WAXY = [none]
UNCX = 1 /HPF
TPO4 = [none]
CAOX = [none]
CAPH = [none]
CACB = [none]
URIC = [none]
LEUC = [none]
CYST = [none]
TYRO = [none]
AMOR = [none]
UNCL = [none]
PC = 1314 /uL